



**CONTAINMENT BOOTH**

**FEATURES AND OPTIONS**  
**CONFIGURATOR**

**CUSTOMER INFORMATION:**

COMPANY: \_\_\_\_\_

CONTACT: \_\_\_\_\_

ADDRESS: \_\_\_\_\_

\_\_\_\_\_

PHONE: \_\_\_\_\_

CELL: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

**ISI REPRESENTATIVE & CONTACT INFO:**

NAME: \_\_\_\_\_

COMPANY: \_\_\_\_\_

PHONE: \_\_\_\_\_

CELL: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

**END LOCATION INFORMATION:**

COMPANY: \_\_\_\_\_

CONTACT: \_\_\_\_\_

ADDRESS: \_\_\_\_\_

\_\_\_\_\_

PHONE: \_\_\_\_\_

CELL: \_\_\_\_\_

E-MAIL: \_\_\_\_\_

**PROJECT TIMEFRAME:**

- ◇ IMMEDIATE NEED (WITHIN 6 MONTHS)
- ◇ CURRENT NEED (6 MONTHS TO ONE YEAR)
- ◇ FUTURE NEED (MORE THAN ONE YEAR)
- ◇ BUDGETARY
- ◇ OTHER, PLEASE SPECIFY: \_\_\_\_\_

**PROJECT INFORMATION:**

- STANDARD FEATURES:**
- ◆ 304L STAINLESS STEEL
  - ◆ PRE-WIRED MODULAR ASSEMBLY
  - ◆ 4" THICK DOUBLE PAN HOLLOW WALL SECTIONS
  - ◆ DAMPERED DISCHARGE GRILLE
  - ◆ FLUSH MOUNTED LED LIGHTING
  - ◆ NON-SPARKING BACKWARD CURVED ALUMINUM IMPELLARS
  - ◆ VARIABLE SPEED CONTROL
  - ◆ CG SCREEN DOWN FLOW AIR DIFFUSERS
  - ◆ MERV-8 EXTRACTION PRE-FILTERS
  - ◆ HIGH VOLUME HEPA EXTRACTION FILTERS
  - ◆ FILTERS, BLOWERS, LIGHTS & CONTROLS ARE ROOM SIDE ACCESSIBLE
  - ◆ VISUAL OPERATIONAL PRESSURE GAUGES
  - ◆ WET LOC. RATED ON / OFF SWITCHES
  - ◆ MUSHROOM HEAD EMERGENCY STOP SWITCH
  - ◆ SINGLE POINT ELECTRICAL CONNECTION



# CONTAINMENT BOOTH FEATURES AND OPTIONS CONFIGURATOR

**URS AVAILABLE?:**

- ◇ YES

**MATERIAL OF CONSTRUCTION:**

- ◇ 304 STAINLESS STEEL #4 FINISH
- ◇ 316L STAINLESS STEEL #4 FINISH
- ◇ SPECIAL FINISH SPECIFY \_\_\_\_\_

**WALL PANELS:**

- ◇ FULL LIGHT CLEAR GLASS WALL PANELS
- ◇ ANTISTATIC CLEAR ACRYLIC WALL PANELS
- ◇ REMOVEABLE CLEAR FLEXIBLE PVC WALLS
- ◇ MODULAR SOFTWALL DESIGNS W/ MODULAR EXTRACTION MODULE

**NEGATIVE PRESSURE OPERATION:**

- ◇ RE-CIRCULATING
- ◇ ONCE THRU / SINGLE PASS
- ◇ FUME EXTRACTION
- ◇ SECONDARY BARRIER SYSTEMS
- ◇ PORTABLE SYSTEMS
- ◇ SMALL BATCH

**STANDARD SIZES:**

**SMALL BATCH:**

- ◇ 24”d x 48”w x 30”h WORK AREA
- ◇ 24”d x 48”w x 36”h WORK AREA
- ◇ 24”d x 60”w x 30”h WORK AREA
- ◇ 24”d x 60”w x 36”h WORK AREA
- ◇ 24”d x 72”w x 30”h WORK AREA
- ◇ 24”d x 72”w x 36”h WORK AREA

**WALK-IN BOOTH:**

- ◇ 4’ x 4’ WORK AREA (78”d x 50”w x 120”h OA)
- ◇ 4’ x 6’ WORK AREA (106”d x 50” w x 120”h OA)
- ◇ 4’ x 8’ WORK AREA (134”d x 50” w x 120”h OA)
- ◇ 6’ x 6’ WORK AREA (106”d x 78”w x 120”h OA)
- ◇ 6’ x 8’ WORK AREA (134”d x 78” w x 120”h OA)
- ◇ 8’ x 8’ WORK AREA (134”d x 102”w x 120”h OA)
- ◇ 8’ x 10’ WORK AREA (162”d x 102”w x 120”h OA)
- ◇ 10’ x 10’ WORK AREA (162”d x 126”w x 120”h OA)

**CUSTOM SIZE OR CONFIGURATION, PLEASE SPECIFY:**

---



---



---



---



---

**OPERATOR EXPOSURE LEVEL:**

- ◇ OEL 2 (100-1,000 µg/m<sup>3</sup>)
- ◇ OEL 3 (10-100 µg/m<sup>3</sup>)

**ISO / FED STD 2019 / EU GMP GRADE:**

- ◇ ISO 8 / CLASS 100,000 / GRADE D
- ◇ ISO 7 / CLASS 10,000 / GRADE C

**ARTICULATING SECONDARY SHIELD:**

- ◇ ACRYLIC
- ◇ POLYCARBONATE
- ◇ XYZ AXIS

**SECONDARY SHIELD GLOVE PORT SIZE:**

- ◇ 8” CIRCULAR
- ◇ 10” CIRCULAR
- ◇ 12” CIRCULAR

**ARTICULATING SECONDARY SHIELD GLOVES:**

- ◇ NATURAL RUBBER
- ◇ NEOPRENE
- ◇ HYPALON (CSM)
- ◇ BUTYL
- ◇ LATEX
- ◇ VITON POLYURETHANE

**DOORS:**

- ◇ SOLID SS DOORS
- ◇ AUTO ROLL UP DOORS
- ◇ STRIP DOORS

**SPECIAL ACCESS DOORS PLEASE SPECIFY:**

---

**FILTRATION:**

- ◇ HEPA 99.99% eff @ 0.3 microns
- ◇ ULPA 99.9995% eff @ 0.12 microns
- ◇ ASHRAE
  - ◇ SPECIFY RATING \_\_\_\_\_
- ◇ CARBON FILTER
  - ◇ SPECIFY RATING \_\_\_\_\_
- ◇ ePTFE HYDROPHOBIC HEPA FILTRATION
- ◇ BAG IN / BAG OUT (BIBO)

**BLOWER SYSTEM OPTIONS:**

- ◇ EXPLOSION PROOF CLASS 1 DIV 1
- ◇ PRESSURE ACTIVATED BLOWER SPEED MODULATION
- ◇ VARIABLE FREQUENCY DRIVE



## CONTAINMENT BOOTH CONFIGURATOR CONTINUED

### INTEGRAL WORK SURFACES:

- ◇ SHELVING \_\_\_\_\_
- ◇ HOOKS \_\_\_\_\_

### INTEGRAL LIGHTING:

- ◇ HIGH INTENSITY BIAX
- ◇ FLUORESCENT
- ◇ INCANDESCENT
- ◇ ULTRA VIOLET
- ◇ VAPOR-TIGHT LED
- ◇ EXPLOSION PROOF CLASS 1 DIV 1 or DIV 2

### ELECTRICAL SERVICES:

- ◇ 120 / 220V 60HZ / 1PH
- ◇ 230 / 460V 60HZ / 3PH

### ELECTRICAL OUTLETS:

- ◇ GFI DUPLEX
- ◇ EXPLOSION PROOF
- ◇ SEALED / STAINLESS STEEL COVER
- ◇ STRIP OUTLET
- ◇ SPECIAL VOLTAGE \_\_\_\_\_
- ◇ SPECIFY AMPERAGE \_\_\_\_\_

### CONTROLS:

- ◇ MAGNEHELIC GAGE
- ◇ MINIHELIC GAGE
- ◇ PHOTOHELIC GAGE
- ◇ DIGIHELIC GAGE
- ◇ PLC & ALARM SYSTEM
- ◇ REMOTE LOCATED CONTROL PANELS
- ◇ AUDIBLE VISUAL ALARM SYSTEM
- ◇ BMS OPERATION

### DRUM INTRODUCTION LIFT:

- ◇ MANUAL HAND CRANK
- ◇ AUTOMATIC ADJUSTMENT
  - ◇ FOOT PEDAL ACTIVATED
  - ◇ PUSH BUTTON ACTIVATED

### ◇ DAMPERED DISCHARGE DUCT CONNECTION

### ◇ EXTRACTION PLENUM CIP SYSTEM

### ◇ CHILLED WATER COIL

### SERVICE PENETRATIONS AND FITTINGS:

- ◇ TRI-CLAMP AND CAP \_\_\_\_\_
- ◇ BULK HEAD FITTINGS \_\_\_\_\_
- ◇ CABLE PENETRATIONS \_\_\_\_\_
- ◇ VALVES & PET COCKS \_\_\_\_\_
- ◇ COMPUTER CONNECTIONS \_\_\_\_\_
- ◇ QUICK CONNECTS \_\_\_\_\_
- ◇ FAUCETS \_\_\_\_\_
- ◇ SPRINKLER PENETRATIONS \_\_\_\_\_
- ◇ OTHER (PLEASE SPECIFY) \_\_\_\_\_

### ◇ FACTORY ACCEPTANCE TESTING (FAT) ◇ SITE ACCEPTANCE TESTING (SAT)

### ◇ EQUIPMENT INSTALLATION

### CLARIFICATION NOTES:

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

---

