

HIGH CONTAINMENT ALPHA/BETA RAPID TRANSFER PORTS



STANDARD FEATURES

- 316 Stainless Steel Construction
- Precision machined components
- Mechanical Interlock
- Alfa Flange
- Alfa Port
- Machined precision Hinge
- Beta Cover
- FDA Approved Silicone Beta Seal
- Beta Flange
- Beta handles
- Beta Container Polypropylene
- No special tools required for container placement or removal.

OPTIONAL FEATURES

- Special Hinges and Doors
- Stainless Steel Beta Containers
- SS Framed Safety Glass Doors
- Autoclave Beta Flange with Containers 316L
- Autoclave Beta Flange with containers Aluminum
- Beta Flange Wrench
- Funnel Assembly
- Bag Ring and Flange Assembly
- Beta flange with Tri Clover Fitting
- Pre Installed service connections

MODEL#	ALFA FLANGE DIA.
RTP-FD-105	105
RTP-FD-190	190
RTP-FD-270	270
RTP-FD-350	350
RTP-FD-350	460

MODEL#	BETA FLANGE DIA.
RTP-BC-105	105
RTP-BC-190	190
RTP-BC-270	270
RTP-BC-350	350
RTP-BC-350	460

CONTAINER LENTHS
300,400,500,600, mm
300,400,500,600, mm
300,400,500,600, mm
300,400,500,600, mm
300,400,500,600, mm

