

# PERSONNEL DECONTAMINATION MISTING SHOWERS



## **STANDARD FEATURES**

- 100% 16 ga. 304 stainless steel construction with #4 finish.
- Fine Atomization Misting nozzles
- Flow Rates 0.023-0.076 gpm
- Direct pressure
- Spray Pattern :Cone shaped fog
- Integral ceiling mounted Wet Location Fluorescent lighting.
- Wet location light switch
- Raised Grated Floor above collection basin
- Manual activated shower cycle.
- 36" CLO doors with full lite safety glass (heavy duty construction)
- Door Closers
- Single point top mounted service connections
- Internally Located removable nozzle panels for easy service
- Bottom Located Drain
- In use indicator light

## **OPTIONAL FEATURES**

- 100% 16 ga. 316L stainless steel construction with #4 finish.
- Polyurethane coated Finish
- Sensor- or push-pad-activated automatic doors.
- Internally located panic circuit for interlock override with audible visual alarm.
- Automatic sensor- activated shower cycle.
- Programmable logic controller for system functions.
- Magnetic Interlock system for doors
- Drenching Shower Head
- Hand held spray ball
- Floor mounted nozzles
- Custom Door configurations
- Drain assist pumping systems
- Custom configurations

## **NOTE:**

- All electrical components are UL listed.
- All electrical wiring is in conformance with NEC & NEMA standards.
- USP-797 compliant
- Isolation Systems Inc. inspection reports available upon request.

MODEL #	OUTSIDE DIM	INSIDE DIM	ELECTRICAL	JETS	DOOR	SHIP WT
MS-44	50"w x 48"d x 102"h	36"w x 48"d x 84"h	120v-60hz-1ph-FLA= 2	45	36" CLO	1300 lb
MS-46	50"w x 72"d x 102"h	36"w x 72"d x 84"h	120v-60hz-1ph-FLA= 2	60	36" CLO	1850 lb
MS-56	60"w x 72"d x 102"h	48"w x 72"d x 84"h	120v-60hz-1ph-FLA= 2	65	36" CLO	1950 lb

All information herein is subject to change without notice  
All shipping weights estimated.

All catalog dimensions are nominal.

