

EQUIPMENT ASSISTANCE JIB LIFT



Shown in Bench Top HEPA Filtered Ventilated Exhausting Enclosure.

STANDARD FEATURES:

- 304 Stainless Steel jib support and cradle construction.
- Push button Lift / Lower.
- 25" Maximum travel.
- Free swing 180° left and right.
- Adjustable cradle straps.
- 100 Lb. capacity.

OPTIONAL FEATURES:

- 316L Stainless Steel jib support and cradle construction.
- Custom sizes available.

MODEL	CRADLE DIMENSIONS
EAJL-105	12"L x 6"W
EAJL-110	18"L x 9"W

**All information herein is subject to change without notice.*

**All catalog dimensions are nominal.*



© 8/2022