

HIGH CONTAINMENT ALPHA / BETA RAPID TRANSFER PORTS (RTP)



STANDARD FEATURES:

- 316 Stainless Steel Construction.
- Precision machined components.
- Mechanical Interlock.
- Alpha Flange.
- Alpha Port.
- Machined precision hinge.
- Beta Cover.
- FDA Approved Silicone Beta Seal.
- Beta Flange.
- Beta handles.
- Beta Container Polypropylene.
- No special tools required for container placement or removal.

OPTIONAL FEATURES:

- Special Hinges and Doors.
- Stainless Steel Beta Containers.
- SS Framed Safety Glass Doors.
- Autoclave Beta Flange with Containers 316L.
- Beta Flange Wrench.
- Funnel Assembly.
- Bag Ring and Flange Assembly.
- Beta flange with Tri Clover Fitting.
- Pre Installed service connections.

MODEL	ALFA FLANGE DIA.
RTP-FD-105	105 mm
RTP-FD-190	190 mm
RTP-FD-270	270 mm
RTP-FD-350	350 mm
RTP-FD-350	460 mm

MODEL	BETA FLANGE DIA.
RTP-BC-105	105 mm
RTP-BC-190	190 mm
RTP-BC-270	270 mm
RTP-BC-350	350 mm
RTP-BC-350	460 mm

BETA CONTAINER LENGTHS
300,400,500,600, mm
300,400,500,600, mm
300,400,500,600, mm
300,400,500,600, mm
300, 400, 500, 600, mm

**All information herein is subject to change without notice.*

**All catalog dimensions are nominal.*



© 7/2021